


Ordering Information




PolyMem Silicone Border Dressings

PolyMem Silicone Border

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
2423 #3 Oval	1.0" x 2.0" (2.5cm X 5.0cm) Pad 2.0" x 3.0" (5.0cm X 7.6cm) Adhesive	50	5 bx / 10 ea	A6212
2435 #5 Oval	2.0" x 3.0" (5.0cm X 7.6cm) Pad 3.5" x 5.0" (8.8cm X 12.7cm) Adhesive	25	5 bx / 5 ea	A6212
2466 #6 Square	3.5" x 3.5" (9.0cm x 9.0cm) Pad 5.2" x 5.2" (13.3cm x 13.3cm) Adhesive	30	2 bx / 15 ea	A6212
2468 #8 Oval	4.0" x 5.7" (10.1cm X 14.6cm) Pad 6.5" x 8.2" (16.5cm X 20.9cm) Adhesive	20	4 bx / 5 ea	A6213
2499 #9 VariSite™	5.0" x 5.1" (12.7cm x 13.1cm) Pad 8.2" x 8.5" (20.9cm x 21.7cm) Adhesive	20	4 bx / 5 ea	A6213
2439 Sacral	4.5" x 4.7" (11.4cm x 12.0cm) Pad 7.2" x 7.8" (18.4cm x 20cm) Adhesive	20	4 bx / 5 ea	A6212

PolyMem® SILICONE BORDER

Our unique silicone border is engineered to stay in place better and more securely than other silicone dressings.



PolyMem Silicone Border Silver

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
1243 #3 Oval	1.0" x 2.0" (2.5cm X 5.0cm) Pad 2.0" x 3.0" (5.0cm X 7.6cm) Adhesive	50	5 bx / 10 ea	A6212
1245 #5 Oval	2.0" x 3.0" (5.0cm X 7.6cm) Pad 3.5" x 5.0" (8.8cm X 12.7cm) Adhesive	25	5 bx / 5 ea	A6212
1248 #8 Oval	4.0" x 5.7" (10.1cm X 14.6cm) Pad 6.5" x 8.2" (16.5cm X 20.9cm) Adhesive	20	4 bx / 5 ea	A6213
1249 Sacral	4.5" x 4.7" (11.4cm x 12.0cm) Pad 7.2" x 7.8" (18.4cm x 20cm) Adhesive	20	4 bx / 5 ea	A6212



PolyMem® and PolyMem Silver® Non-Adhesive Dressings

PolyMem Non-Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5022	1.8" x 1.8" (4.7cm x 4.7cm)	100	5 bx / 20 ea	A6209
5033	3" x 3" (8cm x 8cm)	60	4 bx / 15 ea	A6209
5044	4" x 4" (10cm x 10cm)	60	4 bx / 15 ea	A6209
5055	5" x 5" (13cm x 13cm)	30	2 bx / 15 ea	A6210
5077	6.5" x 7.5" (17cm x 19cm)	30	2 bx / 15 ea	A6210
5124	4" x 12.5" (10cm x 32cm)	12	1 bx / 12 ea	A6210
5244	4" x 24" (10cm x 61cm) Roll	8	2 bx / 4 ea	A6211
5824	8" x 24" (20cm x 60cm) Roll	2	1 bx / 2 ea	A6211


PolyMem Silver Non-Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
1022	1.8" x 1.8" (4.7cm x 4.7cm)	100	5 bx / 20 ea	A6209
1044	4.25" x 4.25" (10.8cm x 10.8cm)	30	2 bx / 15 ea	A6210
1077	6.5" x 7.5" (17cm x 19cm)	15	1 bx / 15 ea	A6211
1124	4.25" x 12.5" (10.8cm x 32cm)	12	1 bx / 12 ea	A6211

PolyMem®

PolyMem® SILVER

The optimal choice when cutting dressings to size and tapeless attachment or roll taping, is your preferred method of application. Use as a primary, secondary or combined primary and secondary dressing.



PolyMem MAX and PolyMem MAX Silver Non-Adhesive Dressings

PolyMem MAX


REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5035	3" x 3" (7.6cm x 7.6cm)	40	4 bx / 10 ea	A6209
5045	4.5" x 4.5" (11cm x 11cm)	20	2 bx / 10 ea	A6210
5088	8" x 8" (20cm x 20cm)	10	2 bx / 5 ea	A6211

PolyMem MAX Silver

1045	4" x 4" (10cm x 10cm)	16	2 bx / 8 ea	A6209
1088	8" x 8" (20cm x 20cm)	10	2 bx / 5 ea	A6211

PolyMem® MAX®

This configuration is a thicker PolyMem pad with high MVTR film for wounds with heavier drainage or when longer wear time is desired. Use as a primary, secondary or combined primary and secondary dressing.



PolyMem WIC, PolyMem WIC Silver, and PolyMem WIC Silver Rope Wound Fillers

PolyMem WIC

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5733	3" x 3" (8cm x 8cm) 4 grams	40	4 bx / 10 ea	A6215
5712	3" x 12" (8cm x 30cm) 16 grams	12	1 bx / 12 ea	A6215
5788	8" x 8" (20cm x 20cm) 28 grams	10	2 bx / 5 ea	A6215

PolyMem WIC Silver


1331	1" x 3" (2.5cm x 7.6cm) 1.3 grams	70	5 bx / 14 ea	A6215
1333	3" x 3" (8cm x 8cm) 4 grams	20	2 bx / 10 ea	A6215
1788	8" x 8" (20cm x 20cm) 28 grams	10	2 bx / 5 ea	A6215

PolyMem WIC Silver Rope

1814	0.4" x 14" (1cm x 35cm) 3 grams	12	2 bx / 6 ea	A6215
------	---------------------------------	----	-------------	-------

PolyMem® WIC®

With no film backing, WIC expands within the wound cavity to fill dead space, keeping the wound bed clean and moisturized while absorbing excess fluid. WIC Rope has a reinforced mesh to ensure one-piece removal.



PolyMem Finger/Toe Dressings

Finger / Toe

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
4401	#1 SM; circ. 1.8" - 2.2" (4.67cm - 5.7cm)	30	5 bx / 6 ea	A6209
4402	#2 MED; circ. 2.2" - 2.6" (5.7cm - 6.72cm)	30	5 bx / 6 ea	A6209
4403	#3 LG; circ. 2.6" - 3.0" (6.72cm - 7.74cm)	30	5 bx / 6 ea	A6209
4404	#4 XL; circ. 3.0" - 3.4" (7.74cm - 8.76cm)	30	5 bx / 6 ea	A6209
4405	#5 XXL; circ. 3.4" - 3.8" (8.76cm - 9.78cm)	30	5 bx / 6 ea	A6210

Finger / Toe Silver

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
1401	#1 SM; circ. 1.8" - 2.2" (4.67cm - 5.7cm)	30	5 bx / 6 ea	A6209
1402	#2 MED; circ. 2.2" - 2.6" (5.7cm - 6.72cm)	30	5 bx / 6 ea	A6209
1403	#3 LG; circ. 2.6" - 3.0" (6.72cm - 7.74cm)	30	5 bx / 6 ea	A6209
1404	#4 XL; circ. 3.0" - 3.4" (7.74cm - 8.76cm)	30	5 bx / 6 ea	A6209
1405	#5 XXL; circ. 3.4" - 3.8" (8.76cm - 9.78cm)	30	5 bx / 6 ea	A6210

This configuration is specially designed for use on fingers and toes. Appropriate for use as a combined primary and secondary dressing or as a secondary dressing.



PolyMem Film Adhesive Dressings

PolyMem Film Adhesive


REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
203	1" x 1" (2.5cm x 2.5cm) Pad 2" x 2" (5cm x 5cm) Adhesive	100	5 bx / 20 ea	A6212
405	2" x 3" (5cm x 7.6cm) Pad 3.5" x 4.5" (8.9cm x 11.4cm) Adhesive	60	4 bx / 15 ea	A6212

PolyMem MAX Film Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
0230	1.25" x 1.5" (3.1cm x 3.8cm) Pad 2.75" x 3" (7cm x 7.6cm) Adhesive	40	4 bx / 10 ea	A6212
606	3.5" x 3.5" (9cm x 9cm) Pad 5.25" x 5.25" (13.3cm x 13.3cm) Adhesive	60	4 bx / 15 ea	A6212
3412	2" x 10" (5cm x 26cm) Pad 3.5" x 11.75" (8.9cm x 29.8cm) Adhesive	12	1 bx / 12 ea	A6213

PolyMem®

The ultra-thin, adhesive-coated urethane film, together with the PolyMem membrane, provides a breathable yet tough, protective, water-resistant dressing of unsurpassed versatility.




PolyMem Surgical

PolyMem Surgical MAX Film Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
0511 Flexible Knee Dressing	3.15" x 9.2" (8cm x 23.35cm) Pad 5.3" x 11" (13cm x 28cm) Adhesive	10	2 bx / 5 ea	A6213
0512 Surgical Dressing	3.5" x 10.5" (9cm x 27cm) Pad 5.5" x 12.75" (14cm x 32cm) Adhesive	10	2 bx / 5 ea	A6213
0481 C-Section Dressing	2.7" x 8.25" (6.85cm x 21cm) Pad 4.7" x 10.2" (12cm x 26cm) Adhesive	10	2 bx / 5 ea	A6213
0230 MAX Film Island Dressing	1.25" x 1.5" (3.1cm x 3.8cm) Pad 2.75" x 3" (7cm x 7.6cm) Adhesive	40	4 bx / 10 ea	A6212
606 MAX Film Island Dressing	3.5" x 3.5" (9cm x 9cm) Pad 5.25" x 5.25" (13.3cm x 13.3cm) Adhesive	60	4 bx / 15 ea	A6212
3412 MAX Film Island Dressing	2" x 10" (5cm x 26cm) Pad 3.5" x 11.75" (8.9cm x 29.8cm) Adhesive	12	1 bx / 12 ea	A6213

PolyMem® SURGICAL

Control inflammatory response for maximum benefit, speed healing, and improve outcomes with a revolutionary dressing for surgical incision sites.



PolyMem Surgical Silver MAX Film Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
1511 Flexible Silver Knee Dressing	3.15" x 9.2" (8cm x 23.35cm) Pad 5.3" x 11" (13cm x 28cm) Adhesive	10	2 bx / 5 ea	A6213
1512 Silver Surgical Dressing	3.5" x 10.5" (9cm x 27cm) Pad 5.5" x 12.75" (14cm x 32cm) Adhesive	10	2 bx / 5 ea	A6213

PolyMem Cloth Adhesive Dressings


PolyMem Cloth Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
7203	1" x 1" (2.5cm x 2.5cm) Pad 2" x 2" (5cm x 5cm) Adhesive	100	5 bx / 20 ea	A6212
7405	2" x 3" (5cm x 8cm) Pad 4" x 5" (10cm x 13cm) Adhesive	60	4 bx / 15 ea	A6212
7606	3.5" x 3.5" (9cm x 9cm) Pad 6" x 6" (15cm x 15cm) Adhesive	60	4 bx / 15 ea	A6212
7031	1" x 1" (2.5cm x 2.5cm) Pad 1" x 3" (2.5cm x 8cm) Adhesive	100	5 bx / 20 ea	A6212
7042	2" x 1.5" (5cm x 4cm) Pad 2" x 4" (5cm x 10cm) Adhesive	100	5 bx / 20 ea	A6212

PolyMem Silver Cloth Adhesive

1766	3.5" x 3.5" (9cm x 9cm) Pad 6" x 6" (15cm x 15cm) Adhesive	30	2 bx / 15 ea	A6212
------	---	----	--------------	-------

A breathable, adhesive-coated cloth backing that provides exceptional comfort for your patient.



Shapes® by PolyMem Oval and Oval Silver Dressings

Shapes by PolyMem Oval Film Adhesive


REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
8023 #3 Oval	1" x 2" (2.5cm x 5cm) Pad 2" x 3" (5cm x 7.6cm) Adhesive	100	5 bx / 20 ea	A6212
8053 #5 Oval	3" x 2" (7.6cm x 5cm) Pad 5" x 3.5" (12.7cm x 8.8cm) Adhesive	30	2 bx / 15 ea	A6212
8086 #8 Oval	4" x 5.7" (10.1cm x 14.6cm) Pad 6.5" x 8.2" (16.5cm x 20.9cm) Adhesive	20	2 bx / 10 ea	A6213

Shapes by PolyMem Oval Silver Film Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
1823 #3 Oval	1" x 2" (2.5cm x 5cm) Pad 2" x 3" (5cm x 7.6cm) Adhesive	100	5 bx / 20 ea	A6212
1853 #5 Oval	3" x 2" (7.6cm x 5cm) Pad 5" x 3.5" (12.7cm x 8.8cm) Adhesive	30	2 bx / 15 ea	A6212
1886 #8 Oval	4" x 5.7" (10.1cm x 14.6cm) Pad 6.5" x 8.2" (16.5cm x 20.9cm) Adhesive	20	2 bx / 10 ea	A6213

by PolyMem®
Shapes

A variety of oval sizes offer a secure fit for wounds on elbows, knees and other contoured areas of the body where the dressing needs to bend and conform.



Shapes® by PolyMem Sacral Dressings


Shapes by PolyMem Sacral Film Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
3709	4.5" x 4.7" (11.4cm x 12.0cm) Pad 7.2" x 7.8" (18.4cm x 20cm) Adhesive	20	2 bx / 10 ea	A6212

Shapes by PolyMem Sacral Silver Film Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
1709	4.5" x 4.7" (11.4cm x 12.0cm) Pad 7.2" x 7.8" (18.4cm x 20cm) Adhesive	20	2 bx / 10 ea	A6212

The heart-shaped dressing is designed specifically for use in the sacral region.



Shapes® by PolyMem Tube Dressings

Shapes by PolyMem Tube Dressings

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5333	2.75" x 2.75" (7cm x 7cm)	30	2 bx / 15 ea	A6209
5335	3.5" x 3.5" (9cm x 9cm)	30	2 bx / 15 ea	A6209



This configuration is pre-cut to enhance patient comfort and help clinicians easily apply the dressing to tube sites. Unlike other dressings used on tube sites, PolyMem dressings do not shed fibers.

Nursicare

Nursicare Sterile Breast Pad

REF	DIMENSIONS	PER CASE	PACKAGING
400	2.5" Diameter (6.30cm Diameter)	60	10 bx / 6 ea

NURSICARE®
by PolyMem®




Nursicare therapeutic breast pads soothe and alleviate pain associated with nursing, while helping to heal cracked or sore nipples.

SportsWrap by PolyMem Dressings

SportWrap by PolyMem Dressings

REF	DIMENSIONS	PER CASE	PACKAGING
9336	3" x 36" (7cm x 91cm) SW Roll	16	2 bx / 8 ea
9536	5" x 36" (12cm x 91cm) SW Roll	8	2 bx / 4 ea
9548	5" x 48" (12cm x 121cm) SW Roll	6	2 bx / 3 ea
9120	2.5" x 12' (6.3cm x 3.6m) SW ST Roll	10	1 bag / 10 ea

by PolyMem
SPORTSWRAP®



Non-sterile dressings created specifically to manage athletic injuries and provide quicker (more effective) management.



www.polymem.com
Made in the USA

Trademarks are owned by or licensed to Ferris. The marks may be registered or pending in several countries. This document is meant for general information. See individual product instructions for use.
© 2019-2022 Ferris Mfg. Corp. MKL-698 R9 1222

FERRIS MFG. CORP.
5133 Northeast Parkway
Fort Worth, TX 76106-1822 USA
Toll Free U.S.A.: 800-765-9636
International: +1 817-900-1301
www.polymem.com